

REDUCING TREATMENT TIME IN SINGLE-TARGET STEREOTACTIC RADIOSURGERY WITH VARIAN IDENTIFY

Nathan Murtha, PhD, MCCPM; Nicolas Ploquin, PhD, FCCPM
Department of Physics & Astronomy, University of Calgary, Canada
Department of Oncology, University of Calgary, Canada
Arthur J.E. Child Comprehensive Cancer Centre, Cancer Care Alberta, Canada



Arthur J.E. Child
Comprehensive
Cancer Centre




UNIVERSITY OF
CALGARY

Introduction


- CBCT between arcs prolongs linac-based SRS treatment (~15-25 mins).
- Increase throughput if we can drop CBCTs safely.
- Single-target SRS suitable for reduced imaging with Identify?
- Identify = stereoscopic surface guidance system.

Single Identify pod observing face




Methods


Pre-arc CBCT, take reference surface



Record IDENTIFY positioning during arc




Post-arc CBCT, compare to Identify means




*18 SRS patients (any met count), 38 arcs

Results

Difference in translations between Identify and CBCT



Difference in rotations between Identify and CBCT



Boxes cover 25th-75th percentile, whiskers 99% of data range. Black line is median, red dot is mean. Dashed red lines are proposed tolerances for Identify.

Discussion

- Commissioning results show sub-mm and sub-degree accuracy (applied shifts, table rotations, camera obstructions, skin tone).
- Patient measurements show good agreement between CBCT and Identify, often within 0.01-0.02 cm and 0.1 deg.
- Proposed Identify tolerances of 0.05 cm and 1.0 deg are suitable for single-target SRS.
- Use of Identify allows dropping of inter-arc CBCTs for single-target SRS, time savings of 15-25 mins per patient.
- Patient and RT experience has improved with shorter treatment times and accurate intra-beam monitoring.

Contact authors at:
nathan.murtha@cancercarealberta.ca

